

CNOT8 antibody

Product Information

Catalog No.:	FNab01812
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


Background

Immunogen information


Immunogen:	CCR4-NOT transcription complex, subunit 8
Synonyms:	CCR4-NOT transcription complex subunit 8 CAF1-like protein (CALIFp) CAF2 CCR4-associated factor 8 Caf1b CNOT8 CALIF POP2
Observed MW:	34 kDa
Uniprot ID :	Q9UFF9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human esophagus using FNab01812(CNOT8 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01812(CNOT8 antibody) at dilution of 1:1000