

CNPY2 antibody

Product Information


Catalog No.:	FNab01814
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8F6
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	canopy 2 homolog(zebrafish)
Synonyms:	Protein canopy homolog 2 MIR-interacting saposin-like protein Putative secreted protein Zsig9 Transmembrane protein 4 CNPY2 MSAP TMEM4 ZSIG9
Observed MW:	21 kDa
Uniprot ID :	Q9Y2B0

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using FNab01814(CNPY2, MSAP antibody) at dilution of 1:500



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab01814(CNPY2, MSAP antibody) at dilution of 1:1000