

CNRIP1 antibody

Product Information

Catalog No.:	FNab01818
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Isoform 1 suppresses cannabinoid receptor CNR1-mediated tonic inhibition of voltage-gated calcium channels. Isoform 2 does not have this effect.

Immunogen information


Immunogen:	cannabinoid receptor interacting protein 1
Synonyms:	CB1 cannabinoid receptor-interacting protein 1 (CRIP-1) CNRIP1 C2orf32
Observed MW:	13-25 kDa
Uniprot ID :	Q96F85

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab01818(CNRIP1 antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01818(CNRIP1 antibody) at dilution of 1:1000