

## CNTFR antibody

### Product Information

Catalog No.:	FNab01820
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds to CNTF. The alpha subunit provides the receptor specificity.

### Immunogen information

Immunogen:	ciliary neurotrophic factor receptor
Synonyms:	Ciliary neurotrophic factor receptor subunit alpha (CNTF receptor subunit alpha, CNTFR-alpha) CNTFR
Observed MW:	41 kDa
Uniprot ID :	P26992

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01820( CNTFR Antibody) at dilution of 1:1000