

CNTN4 antibody

Product Information

Catalog No.:	FNab01822
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Contactins mediate cell surface interactions during nervous system development. Has some neurite outgrowth-promoting activity. May be involved in synaptogenesis.

Immunogen information

Immunogen:	contactin 4
Synonyms:	Contactin-4 Brain-derived immunoglobulin superfamily protein 2 (BIG-2) CNTN4
Observed MW:	113 kDa
Uniprot ID :	Q8I WV2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01822(CNTN4 antibody) at dilution of 1:300