

COG7 antibody

Product Information

Catalog No.:	FNab01831
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for normal Golgi function.

Immunogen information


Immunogen:	component of oligomeric golgi complex 7
Synonyms:	Conserved oligomeric Golgi complex subunit 7 (COG complex subunit 7) Component of oligomeric Golgi complex 7 COG7
Observed MW:	86kd
Uniprot ID :	P83436

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab01831(COG7 antibody) at dilution of 1:100



Y79 cells were subjected to SDS PAGE followed by western blot with FNab01831(COG7 antibody) at dilution of 1:400