

COG8 antibody

Product Information

Catalog No.:	FNab01832
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for normal Golgi function.

Immunogen information

Immunogen:	component of oligomeric golgi complex 8
Synonyms:	Conserved oligomeric Golgi complex subunit 8 (COG complex subunit 8) Component of oligomeric Golgi complex 8 COG8
Observed MW:	68 kDa
Uniprot ID :	Q96MW5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01832(COG8 antibody) at dilution of 1:400