
Collagen Type II antibody

Product Information

Catalog No.:	FNab01837
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	collagen, type II, alpha 1
Synonyms:	Collagen alpha-1(II) chain Alpha-1 type II collagen Collagen alpha-1(II) chain Chondrocalcin COL2A1
Observed MW:	141 kDa
Uniprot ID :	P02458

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung tissue slide using FNab01837(COL2A1 Antibody) at dilution of 1:50



Mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab01837(COL2A1 antibody) at dilution of 1:500