

COL8A1 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab01847 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information


| | |
|--------------|--|
| Immunogen: | collagen, type VIII, alpha 1 |
| Synonyms: | Collagen alpha-1(VIII) chain Endothelial collagen Vastatin COL8A1 C3orf7 |
| Observed MW: | 73 kDa |
| Uniprot ID : | P27658 |

Application

| | |
|-----------------------|--|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IF, IP |
| Recommended dilution: | WB: 1:500-1:3000; IP: 1:500-1:1000; IF: 1:20-1:200 |
| Image: | |



IP Result of anti-Collagen Type VIII
(IP:FNab01847, 4ug; Detection:FNab01847 1:500)
with HepG2 cells lysate 1800ug.



HepG2 cells were subjected to SDS PAGE
followed by western blot with
FNab01847(Collagen Type VIII antibody) at
dilution of 1:500