

COL16A1 antibody

Product Information

Catalog No.:	FNab01848
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	collagen, type XVI, alpha 1
Synonyms:	Collagen alpha-1(XVI) chain COL16A1
Observed MW:	158 kDa
Uniprot ID :	Q07092

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:3000
Image:	

human placenta tissue were subjected to SDS
PAGE followed by western blot with
FNab01848(Collagen Type XVI antibody) at
dilution of 1:1500

