

CFB antibody

Product Information

Catalog No.:	FNab01853
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	complement factor B
Synonyms:	Complement factor B C3/C5 convertase Glycine-rich beta glycoprotein (GBG) PBF2 Properdin factor B Complement factor B Ba fragment Complement factor B Bb fragment CFB BF BFD
Observed MW:	93-100 kDa, 60-69 kDa
Uniprot ID :	P00751

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab01853(CFB Antibody) at dilution of 1:50



human serum tissue were subjected to SDS PAGE followed by western blot with FNab01853(CFB antibody) at dilution of 1:800