

GJB2 antibody

Product Information

Catalog No.:	FNab01857
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	gap junction protein, beta 2, 26kDa
Synonyms:	Gap junction beta-2 protein Connexin-26 (Cx26) GJB2
Observed MW:	26 kDa
Uniprot ID :	P29033

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab01857(Connexin-26 antibody) at dilution of 1:500