

GJA3 antibody

Product Information

Catalog No.:	FNab01862
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	gap junction protein, alpha 3, 46kDa
Synonyms:	Gap junction alpha-3 protein Connexin-46 (Cx46) GJA3
Observed MW:	48 kDa
Uniprot ID :	Q9Y6H8

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:1000
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab01862(Connexin-46 antibody) at dilution of 1:600