

## COQ4 antibody

### Product Information

Catalog No.:	FNab01877
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Component of the coenzyme Q biosynthetic pathway. May play a role in organizing a multi-subunit COQ enzyme complex required for coenzyme Q biosynthesis. Required for steady-state levels of other COQ polypeptides.

### Immunogen information

Immunogen:	coenzyme Q4 homolog
Synonyms:	Ubiquinone biosynthesis protein COQ4 homolog, mitochondrial Coenzyme Q biosynthesis protein 4 homolog COQ4
Observed MW:	32 kDa
Uniprot ID :	Q9Y3A0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab01877(COQ4 antibody) at dilution of 1:2000