

CORO2B antibody

Product Information

Catalog No.:	FNab01884
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in the reorganization of neuronal actin structure.

Immunogen information


Immunogen:	coronin, actin binding protein, 2B
Synonyms:	Coronin-2B Coronin-like protein C (Clipin-C) Protein FC96 CORO2B KIAA0925
Observed MW:	54 kDa
Uniprot ID :	Q9UQ03

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab01884(CORO2B antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01884(CORO2B antibody) at dilution of 1:300