

C1GALT1C1 antibody

Product Information

Catalog No.:	FNab01887
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	C1GALT1-specific chaperone 1
Synonyms:	C1GALT1-specific chaperone 1 C38H2-like protein 1 (C38H2-L1) Core 1 beta1,3-galactosyltransferase 2 (C1Gal-T2, C1GalT2, Core 1 beta3-Gal-T2) Core 1 beta3-galactosyltransferase-specific molecular chaperone C1GALT1C1 COSMC
Observed MW:	36 kDa
Uniprot ID :	Q96EU7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab01887(COSMC antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01887(COSMC antibody) at dilution of 1:800