

COX15 antibody

Product Information

Catalog No.:	FNab01893
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in the biosynthesis of heme A.

Immunogen information


Immunogen:	COX15 homolog, cytochrome c oxidase assembly protein(yeast)
Synonyms:	Cytochrome c oxidase assembly protein COX15 homolog COX15
Observed MW:	44-55 kDa
Uniprot ID :	Q7KZN9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab01893(COX15 antibody) at dilution of 1:100



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab01893(COX15 antibody) at dilution of 1:1000