

COX6A2 antibody

Product Information

Catalog No.:	FNab01903
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.

Immunogen information

Immunogen:	cytochrome c oxidase subunit VIa polypeptide 2
Synonyms:	Cytochrome c oxidase subunit 6A2, mitochondrial Cytochrome c oxidase polypeptide VIa-heart (COXVIAH) Cytochrome c oxidase subunit VIA-muscle (COX VIa-M) COX6A2 COX6A COX6AH
Observed MW:	None
Uniprot ID :	Q02221

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA
Recommended dilution:	None