

## COX6B2 antibody

### Product Information

Catalog No.:	FNab01905
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Connects the two COX monomers into the physiological dimeric form.

### Immunogen information

Immunogen:	cytochrome c oxidase subunit VIb polypeptide 2
Synonyms:	Cytochrome c oxidase subunit 6B2 Cancer/testis antigen 59 (CT59) Cytochrome c oxidase subunit VIb isoform 2 (COX VIb-2) Cytochrome c oxidase subunit VIb, testis-specific isoform COX6B2
Observed MW:	11 kDa
Uniprot ID :	Q6YFQ2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:1000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using FNab01905(COX6B2 Antibody) at dilution of 1:50



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab01905(COX6B2 antibody) at dilution of 1:100