

## COX6C antibody

### Product Information

Catalog No.:	FNab01906
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.

### Immunogen information


Immunogen:	cytochrome c oxidase subunit VIc
Synonyms:	Cytochrome c oxidase subunit 6C Cytochrome c oxidase polypeptide VIc COX6C
Observed MW:	9 kDa
Uniprot ID :	P09669

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab01906(COX6C antibody) at dilution of 1:50



human heart tissue were subjected to SDS PAGE followed by western blot with FNab01906(COX6C antibody) at dilution of 1:400