

## COX7A2L antibody

### Product Information

Catalog No.:	FNab01908
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the regulation of oxidative phosphorylation and energy metabolism(By similarity).  
Necessary for the assembly of mitochondrial respiratory supercomplex(By similarity).

### Immunogen information


Immunogen:	cytochrome c oxidase subunit VIIa polypeptide 2 like
Synonyms:	Cytochrome c oxidase subunit 7A-related protein, mitochondrial COX7a-related protein Cytochrome c oxidase subunit VIIa-related protein EB1 COX7A2L COX7AR COX7RP
Observed MW:	14 kDa
Uniprot ID :	O14548

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab01908(COX7A2L antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01908(COX7A2L Antibody) at dilution of 1:600