

COX7B antibody

Product Information

Catalog No.:	FNab01909
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport. Plays a role in proper central nervous system(CNS) development in vertebrates.

Immunogen information


Immunogen:	cytochrome c oxidase subunit VIIb
Synonyms:	Cytochrome c oxidase subunit 7B, mitochondrial Cytochrome c oxidase polypeptide VIIb COX7B
Observed MW:	9 kDa
Uniprot ID :	P24311

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab01909(COX7B antibody) at dilution of 1:100



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab01909(COX7B antibody) at dilution of 1:400