

CPLX3 antibody

Product Information

Catalog No.:	FNab01923
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Positively regulates a late step in synaptic vesicle exocytosis.

Immunogen information


Immunogen:	complexin 3
Synonyms:	Complexin-3 Complexin III (CPX III) CPLX3
Observed MW:	20-23 kDa
Uniprot ID :	Q8WVH0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-CPLX3 (IP:FNab01923, 3ug;
Detection:FNab01923 1:300) with mouse brain
tissue lysate 6000ug.



mouse brain tissue were subjected to SDS PAGE
followed by western blot with FNab01923(CPLX3
antibody) at dilution of 1:300