

CPNE5 antibody

Product Information

Catalog No.:	FNab01929
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probable calcium-dependent phospholipid-binding protein that may play a role in calcium-mediated intracellular processes(By similarity). Plays a role in dendrite formation by melanocytes(PubMed:23999003).

Immunogen information

Immunogen:	copine V
Synonyms:	Copine-5 Copine V CPNE5 KIAA1599
Observed MW:	68 kDa
Uniprot ID :	Q9HCH3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse embryo tissue were subjected to SDS PAGE followed by western blot with FNab01929(CPNE5 Antibody) at dilution of 1:300