

CPS1 antibody

Product Information

Catalog No.:	FNab01934
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the urea cycle of ureotelic animals where the enzyme plays an important role in removing excess ammonia from the cell.

Immunogen information

Immunogen:	carbamoyl-phosphate synthetase 1, mitochondrial
Synonyms:	Carbamoyl-phosphate synthase [ammonia], mitochondrial Carbamoyl-phosphate synthetase I (CPSase I) CPS1
Observed MW:	150 kDa
Uniprot ID :	P31327

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab01934(CPS1 Antibody) at dilution of 1:1000