

CPXCR1 antibody

Product Information

Catalog No.:	FNab01948
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	CPX chromosome region, candidate 1
Synonyms:	CPX chromosomal region candidate gene 1 protein Cancer/testis antigen 77 (CT77) CPXCR1
Observed MW:	28-30 kDa
Uniprot ID :	Q8N123

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-CPXCR1 (IP:FNab01948, 4ug;
Detection:FNab01948 1:300) with mouse skeletal
muscle tissue lysate 3600ug.



mouse liver tissue were subjected to SDS PAGE
followed by western blot with
FNab01948(CPXCR1 Antibody) at dilution of
1:300