

CPZ antibody

Product Information

Catalog No.:	FNab01949
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Cleaves substrates with C-terminal arginine residues. Probably modulates the Wnt signaling pathway, by cleaving some undefined protein. May play a role in cleavage during prohormone processing.

Immunogen information

Immunogen:	carboxypeptidase Z
Synonyms:	Carboxypeptidase Z (CPZ) CPZ
Observed MW:	None
Uniprot ID :	Q66K79

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab01949 (CPZ antibody) at dilution of 1:100