

CRB3 antibody

Product Information

Catalog No.:	FNab01953
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Involved in the establishment of cell polarity in mammalian epithelial cells. Regulates the morphogenesis of tight junctions.

Immunogen information

Immunogen:	crumbs homolog 3(Drosophila)
Synonyms:	Protein crumbs homolog 3 CRB3
Observed MW:	None
Uniprot ID :	Q9BUF7

Application

Reactivity:	Human
Tested Application:	ELISA
Recommended dilution:	None