

CREB3L1 antibody

Product Information

Catalog No.:	FNab01963
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	cAMP responsive element binding protein 3-like 1
Synonyms:	Cyclic AMP-responsive element-binding protein 3-like protein 1 (cAMP-responsive element-binding protein 3-like protein 1) Old astrocyte specifically-induced substance (OASIS) Processed cyclic AMP-responsive element-binding protein 3-like protein 1 CREB3L1 OASIS
Observed MW:	57 kDa
Uniprot ID :	Q96BA8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab01963(CREB3L1,OASIS antibody) at dilution of 1:50



A549 cells were subjected to SDS PAGE followed by western blot with FNab01963(CREB3L1,OASIS antibody) at dilution of 1:300