

CRIP2 antibody

Product Information

Catalog No.:	FNab01977
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	cysteine rich protein 2
Synonyms:	Cysteine-rich protein 2 (CRP-2) Protein ESP1 CRIP2 CRP2
Observed MW:	23-29 kDa
Uniprot ID :	P52943

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab01977(CRIP2 Antibody) at dilution of 1:1500