

CRISP2 antibody

Product Information

Catalog No.:	FNab01979
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May regulate some ion channels' activity and thereby regulate calcium fluxes during sperm capacitation.

Immunogen information


Immunogen:	cysteine-rich secretory protein 2
Synonyms:	Cysteine-rich secretory protein 2 (CRISP-2) Cancer/testis antigen 36 (CT36) Testis-specific protein TPX-1 CRISP2 GAPDL5 TPX1
Observed MW:	27 kDa
Uniprot ID :	P16562

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab01979(CRISP2 antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab01979(CRISP2 antibody) at dilution of 1:800