

CRLS1-Specific antibody

Product Information

Catalog No.:	FNab01983
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	cardiolipin synthase 1
Synonyms:	Cardiolipin synthase (CMP-forming) (CLS) Protein GCD10 homolog CRLS1 C20orf155 CLS1
Observed MW:	32 kDa
Uniprot ID :	Q9UJA2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab01983(CRLS1-Specific Antibody) at dilution of 1:200