

DPYSL2 antibody

Product Information

Catalog No.:	FNab01986
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	dihydropyrimidinase-like 2
Synonyms:	Dihydropyrimidinase-related protein 2 (DRP-2) Collapsin response mediator protein 2 (CRMP-2) N2A3 Unc-33-like phosphoprotein 2 (ULIP-2) DPYSL2 CRMP2 ULIP2
Observed MW:	62-66 kDa
Uniprot ID :	Q16555

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab01986(DPYSL2 antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01986(CRMP2 antibody) at dilution of 1:500