

## Dpysl3 antibody

### Product Information

Catalog No.:	FNab01988
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	dihydropyrimidinase-like 3
Synonyms:	Dihydropyrimidinase-related protein 3 (DRP-3) Collapsin response mediator protein 4 (CRMP-4) TOAD-64/Ulip/CRMP TUC-4b Dpysl3 Crmp4
Observed MW:	60-70 kDa
Uniprot ID :	Q62952

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab01988(DPYSL3 Antibody) at dilution of 1:50



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab01988(CRMP4 antibody) at dilution of 1:800