

DPYSL5 antibody

Product Information

Catalog No.:	FNab01990
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	dihydropyrimidinase-like 5
Synonyms:	Dihydropyrimidinase-related protein 5 (DRP-5) CRMP3-associated molecule (CRAM) Collapsin response mediator protein 5 (CRMP-5) UNC33-like phosphoprotein 6 (ULIP-6) DPYSL5 CRMP5 ULIP6
Observed MW:	61 kDa
Uniprot ID :	Q9BPU6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab01990(DPYSL5 antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01990(DPYSL5 antibody) at dilution of 1:1000