

CRYBB3 antibody

Product Information

Catalog No.:	FNab02004
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Crystallins are the dominant structural components of the vertebrate eye lens.

Immunogen information

Immunogen:	crystallin, beta B3
Synonyms:	Beta-crystallin B3 Beta-B3 crystallin Beta-crystallin B3, N-terminally processed CRYBB3 CRYB3
Observed MW:	24 kDa
Uniprot ID :	P26998

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse eye tissue were subjected to SDS PAGE followed by western blot with FNab02004(CRYBB3 antibody) at dilution of 1:500