

---

## CSPP1 antibody

### Product Information

Catalog No.:	FNab02026
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in cell-cycle-dependent microtubule organization.

### Immunogen information

Immunogen:	centrosome and spindle pole associated protein 1
Synonyms:	Centrosome and spindle pole-associated protein 1 CSPP1 CSPP
Observed MW:	141 kDa
Uniprot ID :	Q1MSJ5

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab02026(CSPP1 antibody) at dilution of 1:200