

CSTF2T antibody

Product Information

Catalog No.:	FNab02039
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a significant role in AAUAAA-independent mRNA polyadenylation in germ cells. Directly involved in the binding to pre-mRNAs(By similarity).

Immunogen information


Immunogen:	cleavage stimulation factor, 3' pre-RNA, subunit 2, 64kDa, tau variant
Synonyms:	Cleavage stimulation factor subunit 2 tau variant CF-1 64 kDa subunit tau variant Cleavage stimulation factor 64 kDa subunit tau variant (CSTF 64 kDa subunit tau variant) TauCstF-64 CSTF2T KIAA0689
Observed MW:	70 kDa
Uniprot ID :	Q9H0L4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab02039(CSTF2T antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02039(CSTF2T antibody) at dilution of 1:500