

## MIA2 antibody

### Product Information

Catalog No.:	FNab02044
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

Plays a role in the transport of cargos that are too large to fit into COPII-coated vesicles and require specific mechanisms to be incorporated into membrane-bound carriers and exported from the endoplasmic reticulum. This protein is required for collagen VII(COL7A1) secretion by loading COL7A1 into transport carriers and recruiting PREB/SEC12 at the endoplasmic reticulum exit sites.


### Immunogen information

Immunogen:	CTAGE family, member 5
Synonyms:	Melanoma inhibitory activity protein 2 (MIA protein 2) CTAGE family member 5 ER export factor Cutaneous T-cell lymphoma-associated antigen 5 Meningioma-expressed antigen 6/11 MIA2 CTAGE5 MEA11 MEA6 MGEA11 MGEA6
Observed MW:	100-110 kDa
Uniprot ID :	O15320

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000

Image:



human testis tissue were subjected to SDS PAGE followed by western blot with FNab02044(CTAGE5 antibody) at dilution of 1:300