

CTPS1 antibody

Product Information

Catalog No.:	FNab02060
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	CTP synthase
Synonyms:	CTP synthase 1 CTP synthetase 1 UTP--ammonia ligase 1 CTPS1 CTPS
Observed MW:	67 kDa
Uniprot ID :	P17812

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, FC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-CTP synthase (IP:FNab02060, 4ug; Detection:FNab02060 1:500) with HeLa cells lysate 1200ug.



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02060(CTP synthase antibody) at dilution of 1:500