

## C1QTNF1 antibody

### Product Information

Catalog No.:	FNab02064
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	C1q and tumor necrosis factor related protein 1
Synonyms:	Complement C1q tumor necrosis factor-related protein 1 G protein-coupled receptor-interacting protein (GIP) C1QTNF1 CTRP1
Observed MW:	35 kDa
Uniprot ID :	Q9BXJ1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using FNab02064(C1QTNF1 Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab02064(C1QTNF1 antibody) at dilution of 1:500