

C1QTNF4 antibody

Product Information

Catalog No.:	FNab02065
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	C1q and tumor necrosis factor related protein 4
Synonyms:	Complement C1q tumor necrosis factor-related protein 4 C1q/TNF-related protein 4 C1QTNF4 CTRP4
Observed MW:	45 kDa
Uniprot ID :	Q9BXJ3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab02065(C1QTNF4 antibody) at dilution of 1:100



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab02065(C1QTNF4 antibody) at dilution of 1:1500