

CUL1 antibody

Product Information


Catalog No.:	FNab02073
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	cullin 1
Synonyms:	Cullin-1 (CUL-1) CUL1
Observed MW:	90 kDa
Uniprot ID :	Q13616

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab02073(CUL1 antibody) at dilution of 1:50



IP Result of anti-CUL1 (IP:FNab02073, 3ug; Detection:FNab02073 1:300) with HeLa cells lysate 2000ug.



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab02073(CUL1 antibody) at dilution of 1:400