

TXLNG antibody

Product Information

Catalog No.:	FNab02108
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome X open reading frame 15
Synonyms:	Gamma-taxilin Environmental lipopolysaccharide-responding gene protein Factor inhibiting ATF4-mediated transcription (FIAT) Lipopolysaccharide-specific response protein 5 TXLNG CXorf15 ELRG LSR5
Observed MW:	66 kDa
Uniprot ID :	Q9NUQ3

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab02108(CXorf15 Antibody) at dilution of 1:1000