

## CCDC86 antibody

### Product Information

Catalog No.:	FNab02133
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	coiled-coil domain containing 86
Synonyms:	Coiled-coil domain-containing protein 86 Cytokine-induced protein with coiled-coil domain CCDC86 CYCLON
Observed MW:	65 kDa
Uniprot ID :	Q9H6F5

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:1000; IF: 1:20-1:200
Image:	



mouse stomach tissue were subjected to SDS  
PAGE followed by western blot with  
FNab02133(Cyclon/CCDC86 antibody) at dilution  
of 1:500