

PPIB antibody

Product Information

Catalog No.:	FNab02136
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	peptidylprolyl isomerase B(cyclophilin B)
Synonyms:	Peptidyl-prolyl cis-trans isomerase B (PPIase B) CYP-S1 Cyclophilin B Rotamase B S-cyclophilin (SCYLP) PPIB CYPB
Observed MW:	24 kDa
Uniprot ID :	P23284

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab02136(PPIB antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02136(Cyclophilin B antibody) at dilution of 1:1000