

CYFIP2 antibody

Product Information

Catalog No.:	FNab02140
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in T-cell adhesion and p53/TP53-dependent induction of apoptosis. Does not bind RNA.

Immunogen information


Immunogen:	cytoplasmic FMR1 interacting protein 2
Synonyms:	Cytoplasmic FMR1-interacting protein 2 p53-inducible protein 121 CYFIP2 KIAA1168 PIR121
Observed MW:	130 kDa
Uniprot ID :	Q96F07

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab02140(CYFIP2 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02140(CYFIP2 antibody) at dilution of 1:1000