

CYP20A1 antibody

Product Information

Catalog No.:	FNab02150
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	cytochrome P450, family 20, subfamily A, polypeptide 1
Synonyms:	Cytochrome P450 20A1 CYP20A1
Observed MW:	45 kDa, 60 kDa
Uniprot ID :	Q6UW02

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer using FNab02150(CYP20A1 antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02150(CYP20A1 antibody) at dilution of 1:2000