

CYP4Z1-Specific antibody

Product Information

Catalog No.:	FNab02175
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	cytochrome P450, family 4, subfamily Z, polypeptide 1
Synonyms:	Cytochrome P450 4Z1 CYPIVZ1 Laurate 7-monooxygenase CYP4Z1
Observed MW:	59 kDa
Uniprot ID :	Q86W10

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:1000
Image:	



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with FNab02175(CYP4Z1-Specific antibody) at dilution of 1:600