

## PPIC antibody

### Product Information

Catalog No.:	FNab02178
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	peptidylprolyl isomerase C(cyclophilin C)
Synonyms:	Peptidyl-prolyl cis-trans isomerase C (PPIase C) Cyclophilin C Rotamase C PPIC CYPC
Observed MW:	23 kDa
Uniprot ID :	P45877

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab02178(CYPC Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02178(CYPC antibody) at dilution of 1:1000